MS 10310030078
10310030093*

VDRL / SÍFILIS

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ENGLISH

CLINICAL IMPORCLINICAL IMPORTANCE

Syphilis is a human infectious disease caused by a spirochete, the *Treponema pallidum*, transmitted primarily sexually. Other possible routes of transmission are transfusion of infected blood, now virtually eliminated through routine serological screening (congenital syphilis), transmitted *in uterus* treponemas from the infected mother to the developing fetus.

Clinically, after an incubation period ranging from 10 to 90 days, as it is inversely related to the amount of inoculated, occurs, in 85% of patients, the emergence of a chancre, which is a solitary and painless lesion, characterizing primary syphilis. Approximately 4 to 10 weeks after the onset of chancre, symptoms such as weight loss, headache, anorexia, myalgia, arthralgia, malaise, low fever, generalized lymphadenopathy and exanthema (present in 75 to 100% of cases) often appear, which characterizes secondary syphilis. It can last from 1 to 2 years. Without treatment, about one third of patients have tertiary syphilis, which can manifest as gum (15%) cardiovascular syphilis (10%) or neurosyphilis (8 to 10%).

Serological tests for syphilis are classified as non-treponemal, most commonly used for screening, such as VDRL (*Venereal Disease Research Laboratory*) and treponemal tests, used as confirmatory tests for reactive sera in screening tests, such as TPHA (*Pallidum Treponema Hemagglutination*) and FTA-Abs (*Treponema pallidum Hemagglutination*).

The **VDRL** from **WAMA Diagnóstica** is a non-treponemal flocculation test for the diagnosis of syphilis, through the research of antibodies (reagin) in serum, plasma, Cerebrospinal fluid (CSF), with the great advantage over classical VDRL because it consists of a stabilized and ready-to-use suspension.

PRINCIPLE OF THE METHOD

When the antigenic suspension of VDRL is mixed with the sample containing antibodies (reagins), the antigen particles floccule and the result of the reaction is observed under a microscope. The absence of flocculation indicates a negative result.

KIT PRESENTATION

REF55250-F: 250 determinations

1. Antigenic suspension (1x5ml)

2. Instructions for use

REF55250C-F: 250 determinations

1. Antigenic suspension (1x5ml)

2. Positive Control Serum (1x0,5ml)

3. Negative Control Serum (1x1ml)

4. Instructions for use

REF55500-F: 500 determinations

1. Antigenic suspension (2x5ml)

2. Instructions for use

REF55500C-F: 500 determinations

1. Antigenic suspension (2x5ml)

2. Positive Control Serum (1x1ml)

3. Negative Control Serum (1x1ml)

4. Instructions for use

MATERIAL REQUIRED BUT NOT SUPPLIED

Test tubes for dilution and titration

Automatic pipettes

Shelf for pipes and point rack

Containers for disposal of the material

Excavated plate

Saline at 0,9%

PREPARATION AND STABILITY OF REAGENTS

1. ANTIGENIC SUSPENSION: Keep at room temperature before use. Homogenize well before use. Don't freeze.

2. POSITIVE CONTROL SERUM: Human serum, positive for VDRL, diluted in saline phosphate buffer, contains 0,095% sodium azide as a preservative.

3. SERUM NEGATIVE CONTROL: Human serum, negative for VDRL, diluted in saline phosphate buffer, contains 0,095% sodium azide as a preservative.

The kit maintains the same performance after first use and is stable until the expiry date described on the label, provided it is kept at the indicated temperature (2 - 8°C).

TRANSPORT

The kit is not affected by the transport as long as it is delivered to the recipient within a maximum of 96 hours.

SAMPLES

Serum, plasma or Cerebrospinal fluid (CSF), free of hemolysis, lipemia and contamination. If necessary, samples can be stored for a maximum of 4 to 6 weeks in the freezer at -20°C. Do not inactivate samples.

PROCEDURE

A. Qualitative Test

Objective: for screening and elimination of negative samples.

1. Pipet 50µl of sample or controls serums in the cavity of the excavated plate.

2. Pipet 20µl of homogenized antigenic suspension in the same cavities as the samples and control serums. It is not necessary to mix these two components.

3. Shake the plate for 4 minutes at 180rpm.

4. Immediately after 4 minutes, observe under the microscope.

ATTENTION: to avoid pro-zone effect it is suggested that the qualitative test be performed with pure (non diluted) and diluted sample.

READING RESULTS

Negative Reaction: Absence of aggregates. Homogeneous aspect.

Weakly Positive Reaction: Presence of small dispersed aggregates.

Positive Reaction: Presence of medium and large aggregates.

B. Semi-Quantitative Test

1. Dilute the sample in saline solution at 1/2, 1/4, 1/8, 1/16, 1/32, and more if necessary.

2. Pipette 50µl of each dilution into each well of the excavated plate.

3. Pipet 20µl of homogenized antigenic suspension at each dilution. It is not necessary to mix these two components.

4. Shake the placa for 4 minutes at 180rpm.

5. Immediately after 4 minutes, observe under the microscope.

SAMPLE TITLE: It will be the last dilution where the presence of aggregates is still visualized.

QUALITY CONTROL

The clinical laboratory must have a Quality Program to ensure that its procedures are performed in accordance with the principles of Good Clinical Laboratory Practice (GLP).

It is recommended for each battery of tests, to evaluate quality control samples, it is suggested to evaluate at least one reagent sample and one non-reagent.

TEST PERFORMANCE

Sensitivity

In 64 truly positive samples, no false negative results were found, so sensitivity was 100%.

Specificity

In 52 truly negative samples no false positive results were found, so the specificity was 100%.

General Panel		Commercial Kit		
		Positive	Negative	Total
WAMA Kit	Positive	64	0	64
	Negative	0	52	52
	Total	64	52	116
Sensitivity				100%
Specificity				100%

PRECISION

Assay

They were used seven three different tests, with three different batches, to evaluate the reproducibility of the test. Negative and positive values were correctly identified 100% of the time.

Intra-assay

Using the same seven samples were performed three tests under the same conditions with three different lots of VDRL to verify the repeatability of the test. Negative and positive values were correctly identified 100% of the time.

INTERFERING

False-positive reactions with VDRL may occur in conditions such as immunizations, infections, pregnancy, malaria, autoimmune diseases (systemic lupus erytheatosus, etc.), malignant diseases, etc. If VDRL is positive, a confirmation test, specific to treponema, should be performed. WAMA Diagnosticcproduces and distributes FTA-ABs (Immunofluorescence), TPHA (Hemagglutination).

PRECAUTIONS AND WARNINGS

- Only for in vitro diagnostic use.
- Carefully read the instructions for use before performing the test.
- The expiration date corresponds to the last day of the month marked on the kit box label.
- Store the reagents between 2 - 8°C. Do not freeze.
- VDRL antigenic suspension contains thorosal as a preservative that can be toxic if ingested.
- As 0,095% sodium azide is used as a preservative, the disposal of reactives should be accompanied by large volumes of water to avoid accumulation of azide residues in the pipelines, as it can react with lead and copper forming highly explosive salts. In addition, azide is toxic when ingested.
- Sera used in the preparation of controls were tested, with negative results for anti-HCV antibody, hepatitis b surface antigen (HbsAg) and anti-HIV. However, as no diagnostic method offers complete safety, from the absence of these and other infectious agents, it is recommended to treat the control sera as potentially infectious materials.
- Do not use kit components after the expiration date.
- Do not replace components of this kit with that of other manufacturers, nor use different batch components and codes.
- Do not use the reagent when it has a visual characteristic in disagreement with the one specified in the FISPQ of the product.
- Maintain the usual safety care when handling the reagent, all samples should be handled as potentially infectious materials.
- Dispose of all material in accordance with local regulations.
- Use Good Laboratory Practices (LBPs) in the preservation, handling and disposal of materials.

MANAGEMENT RISK

WAMA Diagnóstica, after reviewing and criticaland detailed analysis of all known and/or predicted hazards, concludes that all risks associated with the Products of the MISCELLANEOUS Line were evaluated, that risk reduction measures have been implemented and that the line's products do not present risks greater than the benefits obtained from their use; and that, if used by qualified and trained professionals, aware of the precautions described in the products, will perform their functions with quality, safety and efficacy.

WARRANTY TERM

The **WAMA Diagnóstica** guarantees the exchange of this diagnostic set, provided that it is within the expiration date and is proven by its Technical Advisory That there were no failures in the execution, handling and conservation of this product. The **WAMA** and their distributors are not responsible for failures in kit performance under these conditions.

IMPORTANT WARNING

The **WAMA Diagnóstica** and its distributors are not responsible for any implications arising directly or indirectly from results obtained from the incorrect use of this product. Once tests are performed in an environment outside the control of the manufacturer and distributor, these may be affected by environmental factors and/ or user error.

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