Imuno-Rápido HIV Triline

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ENGLISH

CLINICAL IMPORTANCE

The human immunodeficiency virus (HIV) attacks immune system cells. The acquired immunodeficiency syndrome (AIDS), described in 1981, is characterized as an advanced stage of infection by HIV which the immune system becomes weakened and the patient is more susceptible to another infections.

The HIV transmission occurs by the percutaneous (with sharp instruments), cutaneous in lacerated skin or by mucous membranes contact with biological fluids of seropositive people, the main fluids are: blood, semen, vaginal secretion and breast milk.

The HIV is a virus with RNA genome, of the family *Lentiviridae*. It belongs to the group of non-cytopathic and non-oncogenic retroviruses that need, to multiply, an enzyme called reverse transcriptase, responsible for the transcription of viral RNA in complementary DNA and that can be integrated into the host genome. The virion of HIV is cylindrical, surrounded by a lipid envelope. Its RNA genome contains 10,000 base pairs and is characterized by the presence of two LTR's (Long Terminal Repeats) and nine regions, three of them are structural proteins encoders (gag, pol and env) and the six others are regulators (tat, rev, nef, vpr, vif).

encode proteins and glycoproteins p-24, p-17, gp-120, gp-41 and reverse transcriptase, proteases and integrases enzymes.

There are two types of HIV, called HIV-1 and HIV-2. The genetic homology between them is approximately 40-45%, determining the code and synthesis of several similar gene products. This fact justifies a certain rate of cross-sorreactivity between both, observed in serological tests.

Genomic variants (subtypes) of HIV-1 and HIV-2 have been described. Thus, HIV-1 isolates are classified into two groups: M (major) and O (outlier), with genetic variability of up to 30% in the env segment. In group M, nine subtypes are identified (A, B, C, D, E, F, G, H and I) and in group O only one. In relation to HIV-2, five subtypes are described (A, B, C, D and E).

The **Imuno-Rápido HIV Triline** is an immunochromatographic test that determines qualitatively the presence of antibodies against recombin proteins corresponding HIV-1 gp-41 and gp-120 glycoproteins of HIV-1 and gp-36 glycoprotein of HIV-2 in serum, plasma and whole blood.

PRINCIPLE OF THE METHOD

Anti-HIV-1 and/or anti-HIV-2 antibodies, if present in the sample, bind to recombinant antigens conjAnti-HIV-1 and/or anti-HIV-2 antibodies, if present in the sample, bind to recombinant antigens conjugated with colloidal gold forming an antigen-antibody complex. It flows through the nitrocellulose membrane of the test plate, and it binds to antigens impregnated in test areas 1 and/or 2, determining the appearance of a colored band.

In the absence of anti-HIV-1 and/or anti-HIV-2 antibodies there will be no appearance of the colored band in areas 1 and/or 2. An immobilized control reagent on the test plate membrane will determine the emergence of colored band, its presence demonstrates that the reagents are functioning properly (control area "C").

KIT PRESENTATION

REF 6220103L-R (10 determinations)

- Test device: 10 units
- Diluent: 1 x 3ml
- Instructions for use

REF 6220203L-R (20 determinations)

- Test device: 20 units
- Diluent: 1 x 3ml
- Instructions for use

REF 6220253L-R (25 determinations)

- Test device: 25 units
- Diluent: 1 x 3ml
- Instructions for use

REF 6220403L-R (40 determinations)

- Test device: 40 units
- Diluent: 2 x 3ml
- Instructions for use

REF 6220203LE-R (20 determinations)

- Test device: 20 units
- Diluent: 1 x 3ml
- Collection Kit Accessories: 20 units (disposable capillary/pipette + lancet)
- Instructions for use

REF 6220253LE-R (25 determinations)

- Test device: 25 units
- Diluent: 1 x 3ml
- Collection Kit Accessories: 25 units (disposable capillary/ pipette + lancet)
- Instructions for use

MATERIAL REQUIRED BUT NOT SUPPLIED

Automatic pipette

Disposable tips

Container for disposal ofmaterial.

Timer

PREPARATION AND STABILITY OF REAGENTS

1. TEST PLATE: should be kept at a temperature between 2 - 30°C. **Do not freeze.** Let the strip reach room temperature before performing the tests, if stored in a refrigerator.
2. SAMPLE THINNER: ready for use. Contains 0.095% sodium azide as a preservative. Stable between 2 - 30°C until the expiration date. **Do not freeze.** Let it reach room temperature before performing the tests, if stored in a refrigerator.

Note: The Kit maintains the same performance after first use and it is stable until the expiration date described on the label, if the storage temperature (2 - 30°C) is respected.

TRANSPORTATION

WAMA Diagnóstica declares that temperature variations during products transport do not affect their quality, as long as the 144-hour interval is respected for products of the Imuno-Rápido line.

After the transport time, the products must be stored as indicated on their label (2 to 30°C).

SAMPLES

Use serum, plasma or whole blood free of hemolysis, lipemia and contamination. Use freshly collected whole blood samples in heparin, EDTA or sodium citrate tubes. Check the correct proportion of sample and liquid anticoagulant to avoid excessive dilution of the sample, which may cause false negative results.

Additionally, the whole blood may be collected by finger prick (Check the procedure below).

It is recommended to use fresh samples. If this is not possible, samples should be stored in a refrigerator between 2 - 8°C for up to 48 hours. Do not freeze whole blood samples.

Serum and plasma samples must be used within 72 hours after collection, if stored in a refrigerator at 2 to 8°C. For long periods of storage, samples should be kept in the freezer at –20°C.

Avoid repeated freezing and thawing cycles as this may cause false results.


Diluted sample may cause false-negative result.

SAMPLE COLLECTION BY FINGER STICK

1. Cleaning the fingertip where the blood will be collected using 70% alcohol.
2. Strongly press the lancet tip on the side of the clean finger to pierce it. Then squeeze your finger tightly until it forms a large drop of blood.
3. Touch the capillary/pipette on the formed drop and collect the whole blood. Avoid the formation of air bubbles in this process.

PROCEDURE


- Allow the test plate (1) to reach room temperature before removing it from the laminated envelope.
- Remove the test plate from the laminated envelope immediately before use;
- Place the test plate on a flat, dry surface.
- Pipette the sample (10µL of serum or plasma, or 20µL of whole blood).


5. Dispense 3 drops (100µl) of the diluent solution (2) into the sample well of the test plate
6. Read the results between 20 and 30 minutes. Do not consider results read after that time.

INTERPRETATION OF THE RESULTS


NEGATIVE: Only one colored band will appear in the control area (C).


HIV-1 POSITIVE: Two colored bands will appear, one in the control area (C) and one in test area 1.

HIV-2 POSITIVE: Two colored bands will appear, one in the control area (C) and one in the test area 2.


HIV-1 and HIV-2 INFECTION: Three colored bands will appear, one in the control area (C), one in test area 1 and one in test area 2.

Note: Any color intensity in the test area should be considered a positive result.

INVALID: If no visible band appears in the test (T) and control (C) area, or if no band appears in the control area (C).

NOTE: A blue color band in the control area region (C) will always be present before the test is performed, if it remains after the test and is not replaced by a purple/red color band the test is considered invalid.

Note: The results should be ignored after the given reading time.

3. The expiration date corresponds to the last day of the month marked on the label of the test plate envelope and kit case.
4. It should be avoided to expose the kit to high temperatures as well as directly to the sun.
5. Do not freeze the test plate as this will cause irreversible deterioration.
6. As 0.095% sodium azide is used as a preservative, the disposal of reactivities must be accompanied by large volumes of water to prevent the accumulation of azide residues in the pipelines, as it can react with lead and copper forming highly explosive salts. In addition, azide is toxic if ingested.
7. Allow the reagents to acquire room temperature before start testing.
8. Do not use kit components after the expiration date.
9. Do not replace components of this kit with other manufacturers, or use different batch components and codes.
10. Do not use the reagent when it has a visuallyin disagreement with the one specified in the PRODUCT FISPQ.
11. Maintaining the usual safety care in handling the reagent, all samples should be handled as potentially infectious materials.
12. Dispose of all material in accordance with local regulations.
13. Use Good Laboratory Practices (GLP) in the preservation, handling and disposal of materials.

RISK MANAGEMENT

WAMA Diagnóstica, after reviewing and critical and detailed analysis of all known and/ or predicted hazards, concludes that all risks associated with the Imuno-Rápido Line Products were evaluated, what risk reduction measures were implemented and that the line's products do not present risks greater than the benefits obtained from their use; and that, if used by qualified and trained professionals, aware of the precautions described in the products, will perform their functions with quality, safety and efficacy.

TERMS OF WARRANTY

WAMA Diagnóstica guarantees the exchange of this diagnostic set, provided that it is within the expiration date and is proven by its Technical Advisory that there were no failures in the execution, handling and conservation of this product.

IMPORTANT WARNING

WAMA Diagnóstica and its distributors are not responsible for errors in the results obtained by this product, since the tests were performed in an environment outside the control of the manufacturer and distributor, which may be affected by environmental factors and/or user error.

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BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAFIA

1. CENTER FOR DISEASE CONTROL AND PREVENTION. First 500.000 AIDS cases in United States, 1995. **MMWR**, **44/46**: p.849-853, 1995.
2. DAFFON, A. et al Rapid simple and reliable doctors office test for antibodies to human immunodeficiency virus-1 in serum. **Clin. Chem.**, v. 36: p. 1312-1316, 1990.
3. MARKOVITZ, D.M. Infection with the human immunodeficiency virus type 2. **Ann. Intern. Med.**, v. 118, p. 211-218, 1993.
4. POTTS, K.E. et al. Variabilidade genômica do HIV no Brasil. **AIDS**, v. 7, p. 1191-1197, 1993.
5. PRICE, C.P. et al. Disposable integrated immunoassay devices. In: Price, C.P. and Newman, D.J. editors. **Principles and Practice of Immunoassay**, 2nd ed. Macmillan Reference: p. 581-603, 1997.
6. RATTNER, L. Molecular biology and pathogenesis of HIV infection. Curr. Opinion **Infect**.

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